

## 各確診檢驗單位 TSH 院際品管定量檢驗之中間精密度 (S3)

送檢批次：CHT2012-04-S3、CHT2012-06-S2、CHT2013-04-S3

	Method code	n	Mean	SD	CV
RH01b	01	3	53.2	3.25	6.1%
RH07b	01	3	53.8	1.97	3.7%
RH14 <sup>註1</sup>	01	3	60.5	4.11	6.8%
RH15	01	3	57.5	3.73	6.5%
CL009	01	3	48.4	2.39	4.9%
CL010	01	3	54.2	2.61	4.8%
CL012	01	3	42.4	1.36	3.2%
CL013	01	3	50.7	3.97	7.8%
CL014a	01	3	55.7	3.06	5.5%
CL015	01	3	56.2	3.47	6.2%
Median for Mean, SD and CV			54.0	3.16	5.8%
Range for Mean, SD and CV			42.4 ~ 60.5	1.36 ~ 4.11	3.2% ~ 7.8%
RH01a	03	3	47.2	0.72	1.5%
RH02c <sup>註1</sup>	03	3	43.0	4.17	9.7%
RH06	03	3	47.9	5.77	12.0%
RH07a	03	3	43.0	3.96	9.2%
RH12	03	3	49.6	1.81	3.6%
RH19	03	3	52.2	4.74	9.1%
RH20	03	3	51.5	1.11	2.2%
CL005	03	3	42.4	3.50	8.3%
CL006a	03	3	50.7	0.44	0.9%
CL008	03	3	43.1	1.99	4.6%
CL011 <sup>註2</sup>	03	3	46.7	4.89	10.5%
Median for Mean, SD and CV			47.2	3.50	8.3%
Range for Mean, SD and CV			42.4 ~ 52.2	0.44 ~ 5.77	0.9% ~ 12.0%
Median for All methods			50.7	3.25	6.1%
Range for All methods			42.4 ~ 60.5	0.44 ~ 5.77	0.9% ~ 12.0%

Method code 01 : RIA ( Sandwich-type ) ; Method code 03 : CLIA ( Sandwich-type )

註 1 : 該單位於 CHT2013-04 更換試藥廠牌

註 2 : 該單位自 CHT2012-08 更換試藥廠牌